



NEXEL Co., Ltd.
55 Magokdong-ro, Seoul,
07802 Republic of Korea

+82-2-2088-8886 T
+82-2-2088-8884 F
support@nexel.co.kr E
www.nexel.co.kr W

CERTIFICATE OF ANALYSIS

PRODUCT NAME	Cardiosight®-S
Catalog Number	C-002
Lot Number	C2-6E1


- ▼ The Cardiosight®-S consists of a population of human induced pluripotent stem cell derived cardiomyocytes. NEXEL Co. Ltd. guarantees this product has passed the following test parameters while adhering to the Cardiosight®-S User Guide. Please keep in mind, electrophysiology measures will vary with plating conditions.

PARAMETER	RESULT
Purity / Identity	96.2 % Cardiomyocytes (+cTnT)
Cell Count & Viability	6.8 million cells (64.3%)
Electrophysiology	
Beat rate (BPM)	44.3
FPD (FPDcF)	268 ms (242 ms)
Drug Response	
E4031 (FPDcF increase at 0.01 μ M > 20%)	Pass
Nifedipine (FPDcF decrease at 0.1 μ M > 30%)	Pass
Mexiletine (Spike amplitude decrease at 10 μ M > 50%)	Pass
Sterility	Pass
Mycoplasma	Negative
Bacteria	Negative

The signature below attests that the above tests were performed according to NEXEL SoP.

Quality Control Test by: Jieun An  Date: November 3, 2021

The signature below attests that the production lot has been cleared for sales & distribution.

Approved for Sale by: Dong-Hun Woo  Date: November 3, 2021